

# Obtaining current-voltage characteristic with Keithley 2701 through RS232 on Python

M. R. Soubkovsky

*Master's in Applied Physics and Physics Engineering  
University of Lorraine and CentraleSupélec*

## Abstract

The goal of this study was to measure the current, resistance or voltage with the DAQ System Keithley 2701 through RS232 protocol using Python.

## 1 Device features and commands used

### 1.1 Keithley 2701

### 1.2 Technical Specifications

Keithley 2701 (Fig. 1) is a DAQ System (Data Acquisition System), in other words it is an equipment that allows the user to make specific measurements. The device is a  $6^{1/2}$ -digit high-performance multimeter/data acquisition system. It can measure voltage (DC and AC), current (DC and AC), resistance (2- and 4-wire), temperature (thermocouple, thermistor, and 4-wire RTD), frequency and period, and test continuity.[1]



Figure 1: Keithley 2701. (Source:<https://www.distrelec.de>)

The user is able to communicate with the device through either Ethernet and RS232.

### 1.3 Coding Interface

#### 1.3.1 General Commands

The general commands in the Keithley 2701 follow the IEEE Std 488.2. These common commands that are supported by the SMU are listed in Table 1. Although commands are shown in uppercase, common commands are not case sensitive and either uppercase or lowercase can be used. Note that although these commands are essentially the same as those defined by the IEEE Std 488.2 standard, the Series 2700 does not strictly conform to that standard.[1]

As an example, on Python the code for obtaining the device identification `*IDN?` would look like:

```
print(keithley.query("*IDN?"))
```

#### 1.3.2 Device Specific Commands

The specific commands that were used to operate the equipment are listed in Table 2.

Code	Name	Description
*IDN?	Identification query	Gives the identification tag of the device
*RST	Reset command	Returns the Series 2700B to default conditions
*TST?	Self-test query	Returns a 0
*CLS	Clear status	Clears all event registers and Error Queue
*TRG	Trigger command	Generates the trigger.EVENT.ID trigger event for use with the trigger model.
*OPC	Operation complete command	Set the Operation Complete bit in the Standard Event Register after all pending commands, including overlapped commands, have completed

Table 1: General Commands

## 2 First tests and communication through MAX

In order to assure that the connection with the device is well established, a series of primary tests is performed. For these tests we used the software Measurement & Automation Explorer (MAX) provided by National Instruments (Fig. 2).

MAX Visa Test Panel (Fig. 3) provides the user with a simple GUI to connect with the device and execute the first series of tests and check if the connection is well established. An extra step that needs to be taken on the Keithley 2701 that is not crucial on the 2602B is setting a Termination Character. On MAX the correct Termination Character can be tried on the VISA Test Panel (Fig. 4).

The identification number (Table 1) is usually the first test that was executed to make sure the connection was well made. As example is shown on Figure 5.

## 3 Python Program

In order to explain the whole program, first, a short description of the libraries and the functions will be given.

### 3.1 Libraries used

A list of libraries that were imported in Python and necessary to the execution of the code that follows is shown below and a brief description of the library is given. Libraries:

- **matplotlib.pyplot**: Provides a MATLAB-like plotting framework inside matplotlib
- **numpy**: adds support for large, multi-dimensional arrays and matrices, along with a large collection of high-level mathematical functions to operate on these arrays
- **visa**: enables you to control all kinds of measurement devices independently of the interface (e.g. GPIB, RS232, USB, Ethernet)
- **tkinter**: standard Python interface to the Tk GUI toolkit
- **time**: provides various time-related functions

Code	Description
:SYSTem:BEEPer:STATe 0	deactivate machine beep
:format:elements READ	give only one value (specified on :FUNction)
:FUNction 'VOLTage'	specify Volt measurement
:FUNction 'CURRent'	specify Current measurement
:FUNction 'RESistance'	specify Resistance measurement
READ?	return value measured

Table 2: Keithley 2701 Specific Commands. Source:[1]




Figure 2: MAX Interface

## 3.2 Functions

The following list of functions lists the functions that are part of the program with its parameters and specifications.

- **connexion\_choice(connection)**: enables the user to choose the identifier of the protocol of connection and connects with the device.  
**connection:**     *string*,     contains the identifier of the connection
- **close\_all()**: close any connection
- **reset()**: reset the Keithley 2701 to the default state
- **measurement(nb,caract,delay=0)**: sends the voltage to the device and measures the current

**nb:**             *int*,             number of measurements

**caract:**       *string*,       string that contains the type of measurement (used to display on GUI)

**delay:**        *float*,        delay bewteen measurements(in ms)

**RETURNS:**    *nd.array*     voltage and measure array

- **complete\_measure(nb,caract[,delay=0])**: calls other functions in order to automatically make all the measurements




Figure 3: MAX - VISA Test Panel




Figure 4: Setting a Termination Character

**nb:** *int*, number of measurements

**caract:** *string*, string that contains the type of measurement (used to display on GUI)

**delay:** *float*, delay between measurements (in ms)

**RETURNS:** *list* [input voltage(*float*), measure of current in Amps(*float*)]

### 3.3 Code

First of all, the libraries described beforehand need to be imported

```

1 import visa
2 import numpy as np
3 import matplotlib.pyplot as plt
4 from tkinter import Button, Tk, Frame, Entry, Label, StringVar, Checkbutton, BooleanVar,
  ↳ IntVar, Radiobutton

```




Figure 5: MAX - VISA Test Panel \*IDN? for Keithley 2400

```
5 import time
```

Then we use a series of functions to perform individual tasks. These functions are thoroughly described in the Section 3.2.

```
6 def connexion_choice(connexion):
7     """Permet de choisir la connexion
8     \nAllow to choose the connexion"""
9     global keithley
10    try:
11        rm = visa.ResourceManager() #import visa
12        rm.list_resources() #import visa
13        keithley = rm.open_resource(connexion, send_end=True, read_termination='\r')
14        ↪ #creat device connexion
15        print("connexion succeed")
16    except:
17        print("Connexion error, check the connexion (GPIO,RS232,USB,Ethernet) and it's
18        ↪ number")
19        raise StopIteration ("Erreur de connexion. Verifier la connexion
20        ↪ (GPIO,RS232,USB,Ethernet) et son numeros \
21        ↪ \nConnexion error, check the connexion
22        ↪ (GPIO,RS232,USB,Ethernet) and it's number")
23    #-----
24    def close_all():
25        """Coupe la connexion
26        \nClose the connexion"""
27        try:
28            reset() #reset smu
29            keithley.write(":SYSTem:BEEPPer:STATe 0") #desactive le beep/deactivate sound
30            keithley.close() #ferme la connexion/close the connexion
31            print("Connexion closed")
32        except:
```

```

29         print("Closing error")
30         raise StopIteration ("Erreur de fermeture \
31                               \nClosing error")
32     #-----
33 def reset():
34     """reset smu"""
35     try:
36         keithley.write("*RST\r")
37         print("instrument reset")
38     except:
39         print("Reset error")
40         raise StopIteration ("Erreur de reset \
41                               \nReset error")
42
43     #-----
44 def measurement(nb,caract,delay=0):
45     """Envoi des tensions et mesure des courants. Enregistre les mesures dans la
46     ↪ variable mesure
47     \nSend tensions and measure currents. Save measures in measure variable"""
48     global tension_input,measure
49     measure=[] #creat array for futur measures
50     tension_input=[]
51     i=0
52     keithley.write(":format:elements READ") #give only value
53     if caract=='Um(V)':
54         keithley.write(":FUNCTION 'VOLTage'") #change to volt measurement
55
56     elif caract=='I(V)':
57         keithley.write(":FUNCTION 'CURRent'")
58
59     elif caract=='Resistance(Ohm)':
60         keithley.write(":FUNCTION 'RESistance'")
61
62     print('\nWrite the real voltage or write any string to stop the measurement.\nWrite
63     ↪ nothing in order to keep the theoretical value')
64     for x in tension_input_temp: #loop on tensions values
65         i+=1
66         print("\nSet voltage to %s(V)" % x)
67         try:
68             real_value=float(raw_input('Input voltage: ') or x)
69         except:
70             print("Measure stopped")
71             break
72         tension_input.append(real_value)
73
74         time.sleep(delay) #delay not needed
75
76     #     y=x
77     y=float(keithley.query("READ?")) #read current value I(A)
78
79     #
80     print('Us(V)\t%s\tpoints' % caract) #display tension,current,points,time

```

```

77         print ("%s\t%s\t%s/%s" % (real_value,y,i,nb)) #display
           ↪ tension,current,points,time
78
79         measure.append(y) #add the current value to a array regrouping all measures
80     return tension_input,measure #return the voltage and current array
81     #-----
82 def complete_measure(nb,caract,delay=0):
83     """Fonction principale prennant les valeurs de tkinter en entrant. Trace les mesures
           ↪ et les sauvegardes.
84     Principal fonction taking tkinter value as input. Plot measures and save them."""
85     global tension_input,measure
86
87     tension_input,measure=measurement(nb,caract) #envoi les tensions choisies et mesure
           ↪ les courants associées
88
89     #obsolete ?
90     if len(tension_input) != len(measure):
91         tension_input=tension_input[0:-(len(tension_input)-len(measure))] #reduce input
           ↪ to measure size
92     if measure:
93         plt.figure(num='Measure ') #plot different figures according to a specific name
94         # plt.clf() #clear the graph to avoid superposing data from the same set (can be
           ↪ deactivated if need to superpose)
95         plt.title("Measure ")
96
97         if caract=='Resistance(Ohm)':
98             plt.ylabel(r'Resistance$(\Omega)$')
99         elif caract=='Um(V)':
100             plt.ylabel(r'$U_m(V)$')
101         else:
102             plt.ylabel(caract)
103         plt.xlabel(r'$U_s(V)$')
104         plt.plot(tension_input,measure, '+', label='Measure ') #display
           ↪ current(input_tension) with dots
105         plt.legend() #add legend to the graph (take label from plot)
106         plt.savefig('Measure %s.svg' % (caract), format='svg', dpi=1000,
           ↪ bbox_inches='tight') #save the graph in a vector file
107         plt.show() #plot data
108
109         np.savetxt('Measure %s.csv' %
           ↪ (caract),np.transpose((tension_input,measure)),delimiter="\t") #save data on
           ↪ a binary file

```

---

After defining the functions, we attribute values to the variables, raising an error if the voltage passes the chosen limit:

---

```

98 connexion='COM1'
99 nb=11          #nb de mesures
100 delay=0        #time between applying voltage and measuring current (not needed)
101 smux='temp'

```

```

102     caract='Um(V)'
103
104     connexion_choice(connexion) #try to connect the device using GPIB or RS232
105     reset() #reset smu
106     keithley.write(":SYSTem:BEEPer:STATe 0") #deactivate beep

```

---

And finally we execute the TKinter code:

---

```

114     def compute():
115         """Fonction utilisée par tkinter pour commander l'instrument
116         \nFonction use by tkinter to pilote the instrument"""
117         msg_floats["text"] = ""
118         msg_nb["text"] = ""
119         msg_high_volt["text"] = ""
120         msg_measure["text"] = ""
121         msg_nb["text"] = ""
122
123         smux=str(smux_entry.get()) #return the smux value in the tkinter entry
124         if measure_choice.get()==0: #see if can measure resistance and volt or resistance
125             ↪ and current                                #if so, change into "if volt and current true"
126             ↪ then change loop measure and caract
127             caract='Um(V)'
128         elif measure_choice.get()==1:
129             caract='I(V)'
130         elif measure_choice.get()==2:
131             caract='Resistance(Ohm)'
132
133         try:
134             volts_min=float(volt_min_entry.get()) # min voltage
135             volts_max=float(volt_max_entry.get()) # max voltage
136             if abs(volts_max)>10 or abs(volts_min)>10:
137                 txt_high_volt="abs(volt) <=10"
138                 msg_high_volt["text"] = txt_high_volt
139                 print(txt_high_volt)
140             else:
141                 nb=int(point_number_entry.get()) #measures nb
142                 if nb<1:
143                     txt_nb="nb>0"
144                     msg_nb["text"] = txt_nb
145                     print(txt_nb)
146                 else:
147                     #             delay=float(delay_entry.get()) #not needed
148                     #             print(smux,volts_min,volts_max,nb,delay) #used to debug
149                     try: #issues with try: hide internals errors
150                         complete_measure(volts_min,volts_max,nb,caract,smux=smux)
151                         txt_measure="Measures done"
152                         msg_measure["text"] = txt_measure
153                         print(txt_measure)
154                     except:
155                         txt_measure="Error from measurement detected"
156                         msg_measure["text"] = txt_measure

```



```

154         print(txt_measure)
155         reset() #cause False positive floats error if using without
            ↪ instrument and previous commands disable
156     except:
157         txt_floats="floats"
158         msg_floats["text"] = txt_floats
159         print(txt_floats)
160
161     root = Tk() #used to creat user interface
162     frame = Frame(root)
163     root.title("KMT - Keithley Measurement Tool") #[Marcel]: I changed the labels just a
        ↪ bit
164     frame.pack()
165
166     L0 = Label(frame, text="Measure's name:") #fixed text
167     L0.grid(row=0, column=0)
168     smux_entry = Entry(frame, textvariable=StringVar(frame, value=smux), bd =2, width=7)
        ↪ #stringvar is used to have default values
169     smux_entry.grid(row=0, column=1) #grid is used to position items on the interface
170
171
172
173     msg_floats = Label(frame, text="")
174     msg_floats.grid(row=2, column=2)
175
176     L3 = Label(frame, text="Number of points")
177     L3.grid(row=3, column=0)
178     point_number_entry = Entry(frame, textvariable=StringVar(frame, value=nb), bd =2,
        ↪ width=7)
179     point_number_entry.grid(row=3, column=1)
180
181     msg_nb = Label(frame, text="")
182     msg_nb.grid(row=3, column=2)
183
184     measure_choice= IntVar()
185     measure_choice.set(1)
186     options = [
187         ("Voltage(V)",
188         ("Current(A)",
189         ("Resistance()")
190     ]
191     Check_button=Label(frame, text="Choose measurement:")
192     Check_button.grid(row=4, column=0)
193
194     for val, language in enumerate(options):
195         option_button=Radiobutton(frame,
196             text=language,
197             variable=measure_choice,
198             value=val)
199         option_button.grid(row=val+4, column=1)
200

```




Figure 6: GUI TKinter with measurement options

```

201 compute_button = Button(frame, text="Measure", width=14, command=compute) #button used
    ↪ to get all values and start measures
202 compute_button.grid(row=8, column=1)
203
204 msg_measure = Label(frame, text="")
205 msg_measure.grid(row=8, column=2)
206
207 root.iconbitmap(default='crystal_oscillator1600.ico')
208
209 root.mainloop() #instance looping until closed
210 close_all() #reset and close connexion with the instrument

```

The GUI can be seen on Fig. 6.

After the clicking on "Compute", a csv file with the data for the voltage and current, and a pdf chart are saved on the same folder as the python code; and the chart is shown with matplotlib.pyplot.

## 4 Practical Tests

A test was performed to obtain the current-voltage characteristic of a PN diode. The result can be seen on Figure 7.




Figure 7: Practical test performed to obtain the Current-Voltage characteristic for a PN diode with the program

## 5 GitHub Depository

This project is public under MIT License on GitHub. The depository is accessible through this link <https://github.com/marcelrsoub/keithley-visa-measurements>.

## References

[1] Series 2700 System SourceMeter Manual.